

ABI-500CL

Advanced – Multi-level

simpleABI[™]
with CUFF:LINK[™]



ABI-500CL Highlights

- ✦ **PC-based system**
 - ✦ Easy staff training - simple to learn and use
 - ✦ Point & Click reporting right on the computer
- ✦ **Cuff-Link**
 - ✦ 8MHz Doppler
 - ✦ Automated cuff selector – no moving tubes
 - ✦ Automated inflation & deflation
- ✦ **8.5x11 Standard Reports**
 - ✦ No taping or pasting waveform strips
- ✦ **Data files easily attached to EMR records**
- ✦ **2-year warranty**

Tests Performed

- ✦ simplePVR[™] Exam
- ✦ ABI Exam
- ✦ ABI with Toe (TBI) Exam
- ✦ 3-cuff Segmental Exam
- ✦ 3-cuff with Toe Segmental Exam
- ✦ 4-cuff Segmental Exam
- ✦ 4-cuff with Toe Segmental Exam

ABI-500CL Kit – what's included:

- ✦ simpleABI Reporting Software
- ✦ 8 port Cuff-Link[™] - cuff selector with remote
- ✦ Sensitive 8MHz Doppler
- ✦ PPG Probe
- ✦ 10 cuffs (2x10cm, 4x12cm, 2x12cm long, 2x2.5cm)

ABI-500CL Options

- ✦ Computer
- ✦ DICOM Plug-in
- ✦ Mobile Carry Bag
- ✦ Additional cuff sizes

Newman  Medical

Employee Owned & Operated

www.newman-medical.com/abi-500cl/

www.newman-medical.com/information/

800-267-5549

WWW.NEWMAN-MEDICAL.COM