

DigiDop



MADE IN



THE USA



BREAKTHROUGH

SENSITIVITY

XD-EXTENDED DEPTH PROBES

XD
eXtended Depth
Probes

56% more sensitive
where you need it.

An Honest

5 YEAR WARRANTY

Our 5 year warranty covers ALL labor and ALL parts Including probes and cords!




ASK YOUR SALES REPRESENTATIVE TODAY!
FOR QUESTIONS OR SALES SUPPORT CALL

800-267-5549

DD-990R Table Top

Room Filling Sound from the large forward-facing speaker of the DD-990R. Its bright large display is easily seen across the room.



Large, Bright Display: Indicates heart rate, battery status, volume setting and signal strength.

Smart  Recharge System

Our **Premium Rechargeable Batteries** stay charged when not in use and can be recharged more than 2000 times.

A **significant** cost saving and reduction of waste

PROBE VERSATILITY: ALL PROBES CAN BE USED ON ALL OUR DOPPLERS

OB: Our XD probes are particularly sensitive to shallow early fetal hearts without sacrificing sensitivity at depth. Both 2 and 3 MHz are available in waterproof versions as well as standard.

Vascular: In addition to our sensitive and easy to hold 5 and 8 MHz probes, a unique audio PPG probe is available for detecting digital or penile flow.

MAIN UNITS	DESCRIPTION
DD-990R	DigiDop 900R Table Top Doppler
DD-300	Non-Display Digital Doppler
DD-301	Non-Display Digital Doppler with Recharger
DD-700	Display Digital Doppler
DD-701	Display Digital Doppler with Recharger

PROBES	DESCRIPTION
D3	General purpose 3MHz probe - best for early fetal hearts
D2	2MHz probe designed for deeper fetal hearts
D3W	3MHz obstetrical for water applications
D2W	2MHz obstetrical for water applications
D8	General purpose 8MHz probe - peripheral arteries
D5	5MHz probe designed for deeper vessels
DPPG	Audio PPG for TBI and other studies

HOW TO ORDER

For product ordering code, add the probe number to the end of the main unit number with a hyphen.

For example, to order a Display Doppler (DD-700) with 8MHz probe (D8) the order code is DD-700-D8.