



DX12 WIRELESS



DP 12 WIRED

SE-1515

PC ECG Workstation

SE-1515 PC-based ECG

- Powerful software with comprehensive functions
- 9/12/15/16/18-lead ECG configurable
- DE15/DE18 sampling box with convenient sampling button and signal quality indicator
- Multi-user Permission Configuration
- Enhanced powerful technical specifications provides reliable signal quality
- Lead reversal function to easily rectify wrong placement of electrodes
- Auto arrhythmia detection highlights the arrhythmia part on ECG traces automatically
- Event marker enables you to take notes anytime during sampling
- Feature description offers you more confidence for the diagnosis conclusion
- Reports comparison provides you an easy way of comparing ECG traces

- Stress test function (Optional)
- Wireless sampling box upgrade available as SE-1515 DX

Data Management System

- Support all ECGs except SE-1/100
- Bi-directional communication with ECGs
- Multi-user Permission Configuration
- Worklist support
- Patient information and ECG data management
- ECG reports printing on A4 paper
- Massive data sharing and management solution
- Multi-format reports transferrable via Email
- Support multi-format reports: XML/SCP/DICOM/PDF/JPG/BMP/TIF/PNG