

VISION M

- Automatic APR setting of the X-ray head
- Control generator settings from laptop
- Modern, ergonomic design
- Tech-friendly retractable stability leg
- Protective, fixed laptop mount
- Secure DR panel storage bin
- Superior analytics - detailed history reports
- Lightweight, powdercoated aluminum frame
- Dual wraps for power and tether cable management

Smart.

Simple.

Intelligent Workflow.



Your partner in mobile DR solutions.

800.349.5980 | 585.272.1690 | sales@fsimed.com

VISION M

The first portable X-ray cart designed around digital imaging

AVANSE Acquire Software

- Advanced anatomy driven image processing
- Automated image stitching
- Repeat / Reject Analysis
- Simple three step exam workflow
- Complete DICOM integration:
 - Multiple PACS, MWL, CD, Print Patient registration
 - Interactive examination worklist
- Integrated image acquisition
- Anatomy dependent image processing (APR)
- Advanced post-processing and image manipulation
- Image annotation and measurements
- DICOM MWL Examination Search/Load, DICOM MPPS Support
- DICOM image formatting
- Local data storage, DICOM image compression

Digital Flat Panel

Panel: VAREX (Varian)
Platform: Wireless
Scintillator: CsI
Active Area: 14" x 17"
Pixel Array: 2476 x 3072
A/D Conversion: 16 Bits
Nyquist: 139µm
Weight: 7.3 lbs

Mobile DR Power Package

- 2.4 kW Generator
- AC Line 120V
- X-ray Range 40 – 100 KvP
- MAS Range 0.4 – 100/25
- MA Range 16 – 35/5
- Unlimited APR's
- Tube focal spot 1.2 x 1.2 mm



Your partner in mobile DR solutions.

800.349.5980 | 585.272.1690 | sales@fsimed.com